

# TECHNICAL SPECIFICATION SHEET

---

## ALG INDUSTRIAL Foldable Respirator Patriot F-01

### KEY FEATURES:

- 99% BFE (Bacterial Filtration Efficiency) according to ASTM F2101
- 99.7% Filter Efficiency according to ICS Laboratories Inc.
- Lightweight construction promotes longer wear time
- Moldable nose clip will help provide custom fit and secure seal
- One size fits most
- Latex free: Yes
- Model: Foldable



### COUNTRY OF ORIGIN:

Made in the USA with  
Globally Sourced Materials

### MATERIAL COMPOSITION:

- **Outer Layer:** Non-woven fabric
- **Filter:** Electrostatic hot air cotton
- **Lining:** Melt blown non-woven fabric
- **Inner Layer:** Spun bound cloth/non-woven fabric
- **Nose Clip:** Aluminum moldable
- **Nose Foam:** Spandex filament sponge bar
- **Headband:** Nylon string

### DO NOT USE:

- For gases and Vapors, including oil aerosols, etc.
- If modified or used for its unintended purpose
- Expiration date of 5 years after manufacture date

### DO USE FOR:

- Medical clinics and operating rooms
- Laboratory procedures
- Infection control requirements
- Emergency planning and preparedness

# FOLDABLE RESPIRATOR

## PACKAGE INFORMATION

Description	PCS/Box	Box/Carton	Pcs/Pkg	Carton/Pallet	Total Quan
PATRIOT F-01	40	20	800	24	19,200

## DIMENSIONS

