

INDUSTRIAL / NON-SURGICAL

“FOLDABLE” N95 RESPIRATOR

NIOSH Industrial Approved – TC-84A-9275

Part Number: PT-N95F-01

Description: Foldable N95 rated respirator with moldable nose clip

Location: Manufactured in Ohio, USA

Brand: Patriot Mask



KEY FEATURES:

- Meets **NIOSH** inhalation procedure testing TEB-APR-STP-0007.
- ≥95% BFE (Bacterial Filtration Efficiency) according to ASTM F2101.
- ≥95% Filter Efficiency according to ICS Laboratories Inc.
- Lightweight construction promotes longer wear time.
- Moldable nose clip will help provide custom fit and secure seal.
- One size fits most.
- Latex free: Yes
- Model: Foldable.

DO NOT USE:

- For gases and Vapors, including oil aerosols, etc.
- For repeated use following NIOSH guidelines
- If modified or used for its unintended purpose
- > Expiration date of 3 years after manufacture date

MATERIAL COMPOSITION:

- **Outer Layer:** Spun bound cloth/non-woven fabric.
- **Filter:** Electrostatic hot air cotton.
- **Lining:** Melt blown non-woven fabric.
- **Inner Layer:** Spun bound cloth/non-woven fabric.
- **Nose Clip:** Aluminum moldable.
- **Nose Foam:** Spandex filament sponge bar.
- **Headband:** Nylon string.
- **Breathing valve:** N/A

DO USE FOR:

- Medical clinics and operating rooms
- Laboratory procedures
- Infection control requirements
- Emergency planning and preparedness

INDUSTRIAL / NON-SURGICAL

“FOLDABLE” N95 RESPIRATOR

NIOSH Industrial Approved - TC-84A-9275

Part Number: PT-N95F-01
Description: Foldable N95 rated respirator with moldable nose clip
Location: Manufactured in Ohio, USA
Brand: Patriot Mask

PACKAGE INFORMATION

Description	PCS/Box	Box/Carton	Pcs/Pkg	Carton/Pallet	Total Quan
PT-N95F-01	40	20	800	32	25,600

DIMENSIONS

