

INDUSTRIAL / NON-SURGICAL

“MOLDED” N95 RESPIRATOR

Manufactured in USA with globally sourced material.

Part Number: Patriot C-02/CS-02
Description: Molded N95 rated respirator with nylon head bands

Location: Manufactured in Ohio, USA

Brand: Patriot Mask



Made in the USA with globally sourced materials!

KEY FEATURES:

- ≥95% BFE (Bacterial Filtration Efficiency) according to ASTM F2101
- ≥95% Filter Efficiency according to ICS Laboratories Inc.
- Pre-formed for secure fit.
- Nose foam
- Moldable nose clip will help provide custom fit and secure seal.
- Stapled headbands
- Available in Small and Regular
- Latex free: Yes
- Model: Molded Cup

MATERIAL COMPOSITION:

- Outer Layer: Spun bound cloth/non-woven fabric.
- Straps - Polyisoprene
- Staples - Ultrasonically welded
- Nose Clip - Aluminum
- Filter - Polypropylene
- Shell - Polyester
- Inner Layer: Needle punch cotton
- Nose Clip: 1 mm Thick Metal
- Nose Foam: Spandex filament sponge bar
- Headband: Nylon string
- Breathing valve: N/A

DO NOT USE:

- For atmospheres less than 19.5% oxygen
- For gases and Vapors, including oil aerosols, etc.
- If modified or used for its unintended purpose
- Expiration date of 3 years after manufacture date

DO USE FOR:

- Medical clinics and operating rooms
- Laboratory procedures
- Infection control requirements
- Emergency planning and preparedness

INDUSTRIAL / NON-SURGICAL

“MOLDED” N95 RESPIRATOR

Manufactured in USA with globally sourced material.


Part Number: Patriot C-02/CS-02
Description: Molded N95 rated respirator with nylon head bands
Location: Manufactured in Ohio, USA
Brand: Patriot Mask

PACKAGE INFORMATION

Description	PCS/Box	Box/Carton	Pcs/Pkg	Carton/Pallet	Total Quan
Patriot C-02	25	20	500	12	12,000
Patriot CS-02	25	20	500	12	12,000

DIMENSIONS

REGULAR



SMALL

